

Training / Prescreening

Pre-arrival:

1. turn on computer
2. check battery of KAT VR sensors
3. replace clean face cover for HTC Vive
4. clean and oil the treadmill surface
5. leave consent form on the desk

Arrival:

1. sign consent form and payment form (if paid)
2. introduce to study
3. show participants how to walk on the treadmill
4. put on shoe covers and sensors
5. adjust the height of treadmill
6. fasten safety belts around waist and legs
7. inform participants to hold on to safety bar when entering

Start the experiment:

1. turn on the waist sensors and shoe sensors
2. click the KAT VR industry platform on the computer desktop
3. rotate the participant facing the wall direction
4. in the KAT VR platform, click "calibration" and "stationary"
5. click Unity or built .exe files to run experiments

Post-experiment:

1. inform participants to hold on to safety bar when exiting
2. remove sensors and shoe covers
3. turn off waist sensors and shoe sensors
4. turn off computer
5. use alcohol wipe to clean the treadmill surface, desks